

# NEWMED

SINCE  
1993

# ECHOSON

ULTRASOUND SCANNERS



## ALBIT DUPLEX Doppler

**MULTI-APPLICATION  
ULTRASOUND SCANNER**



# ALBIT

ULTRASOUND  
WHICH WILL  
AMAZE YOU...




ECHOSON

**ALBIT** - easy-to-use, applicative and operative ultrasound scanner for general diagnostic applications such as GP, OB-Gyn, Urology, Pediatric and specific diagnostic applications such as Colorectal - ERUS, ~~Private Bi-Plane TRUS~~, Anesthesia, Vascular Surgery, Musculoskeletal and Orthopedics.

## PARTICULAR FEATURES:

- **User friendly multi-mode interface** via touch-screen menu or/and operative keyboard, intuitive, quick and easy operation, full HD 15,6" monitor
- **Remote operation controller** enables to operate device from a distance and control basic and most often used functions (programmable)
- **Duplex Color Doppler (CF)**, Power Doppler (PF) and Pulsed Wave (PW) functions as the extension and support for classic 2D mode
- **Multi-Modus** Measurements performed by touching particular points on image, classic - by track-ball and grid estimation.
- **Modular ergonomic** arrangements trolley (bedside), portable, on-table, wall or ceiling month
- **Operational touch-pen, stylus** can be used to control scanner, even with "dirty" hands /covered with gel or blood/, hygienic separation



## Modular system - full set version



## Modular system - on-table version

# MULTI-FREQUENCY PROBES



linear:  
LA-510 / L40 (5,0~17,0 MHz)



convex:  
CA-255 / R60 (2,5~5,0 MHz)



convex:  
CA-409 / R13 (4,0~10,0 MHz)





micro-convex endovaginal:  
CV-580 / R13 (5,0~10,0 MHz)




anorectal 360°:  
R-510 / R8 (5,0~12,0 MHz)

# APPLICATION BAZAAR



**NEWMED**  
MEDITUONTI OY  
045 157 7336  
newmed@newmed.fi  
www.newmed.fi



**ECHOSON**

+ 48 81 886 36 13  
export@echoson.eu  
www.echoson.eu

*ECHO-SON reserves the right to change design and specification without notice. All rights reserved.*



European Union  
European Regional  
Development Fund

